



# SAM Pelvic Sling II





# SAM Pelvic Sling II

- **Step 1:**
  - Empty the casualty's pockets and remove items attached to his belt.





# SAM Pelvic Sling II

- **Step 2:**
  - Pass the binder behind the thighs and slide it upward to the correct position at the level of the greater trochanters.





# SAM Pelvic Sling II

- **Step 3:**
  - **Pass the black strap through the buckle and pull to tighten using counter-traction until the device is heard or felt to click.**





# SAM Pelvic Sling II

- **Step 4:**
  - **Maintain tension and attach the Velcro to itself in front.**





# SAM Pelvic Sling II

- **Step 5:**
  - **Secure the ankles - not too tightly. (Note the use of a tourniquet here for this purpose.)**

